



The CT scan and viewer are in a folder on the Desktop. Double click this folder “RuthCTfoot2020” to open it.



That will open the window on the next page.


Double click PBCVIEW.EXE.

It can take a few seconds to launch, so be patient.



A large black and grey window will come up with your case name pre-loaded.

Click or double click (I don't know which) on your name.



Then you will have to wait for the study to come up. There is a blue progress bar on the bottom which will give you an idea of how long it will take... probably 30 to 60 seconds.

Then use the program. Using the trackpad to navigate can be a bit tricky... Call me if you can't figure it out. May need to hold one finger while scrolling with the other to drag and drop.

The program may need a few moments to generate high-resolution images. You may need to be patient...